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WALL Family

NE41 12005-00 v2.0

Information in this user manual is subject to change without notice.

NEAT Electronics AB reserves the right to change or improve their products and to make changes to the content without obligation to notify any person or organization of such changes or improvements.

NEAT Electronics AB is not responsible for any loss of data, income or any consequential damage whatsoever caused.

For more information, details and descriptions, visit our web site:

www.neat-group.com

Declaration of Conformity

Hereby NEAT Electronics AB declares that the radio equipment WALL Family is in compliance with Directive 2014/53/EU. The full text of the EU declaration of conformity is available at the following internet address:

www.neat-group.com/downloads/documentation

© 2012 Legrand

All rights reserved.

Document number: NE41 12005-00 v2.0

Revision date: 2020-04-22

NEAT Electronics AB

Varuvägen 2

SE-246 42 Löddeköping

Sweden

Phone: +46 (0)46 707065

Fax: +46 (0)46 707087

www.neat-group.com

infosweden@neat-group.com

US Notes

FCC ID: 2AGLF1123001

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

1 Important

1.1 Safety notes

- Read instructions prior to use.
- Always test the system per instructions prior to use.
- The product may not be suitable for all persons.
- Check device regularly and replace when necessary.
- Always check the function of the product after making adjustments.
- Our units are NOT intended for any life support device, thus intending a device whose malfunction may result in damage to a life.

1.2 Use

- Use only original parts.
- Do not expose to direct sunlight.
- Keep away from dust, moist and dirt.
- Do not drop, knock, twist or shake the device.
- Do not warm up the device or use it near fire.
- Products in the WALL Family may not be painted.
- For repairs, contact a NEAT dealer.

1.3 Cleaning

- Clean the device with a soft cloth, dampened slightly with mild soapy water.
- Do not clean the device with harsh chemicals, solvents or other corrosive substances.

1.4 Recycling

Dispose of properly. The worn out product must be returned to a recycling facility for proper disposal or returned to NEAT.

2 About WALL Family products

The WALL Family is divided into 11 products in two categories, WIOR and ROOM. Each product with a different hardware setup to accommodate special needs and requests. The terms WALL, WALL units and WALL family all refer to a common denomination for all 11 products.

The products are equipped with different connectors for different functions:

- A key cylinder to lock/unlock the unit from operation.
- Pull cord to trig alarms by pulling the attached string.
- Push button(s) to trig alarms.
- 6.35 mm tele jack for e.g. a pear button.

Some WALL products are also available as "IP44 Ready". This means that combined with an IP44 mounting frame, the unit can be installed in areas where there is a risk of exposure to humidity and splashes according to an IP44 rating.

3 Product variant denominations

Product Button Connector IP44 Ready

WIOR Category products

WIOR	-	-	Yes
KEY	-	Key cylinder	
PULL	-	Pull cord	Yes
PEAR	-	6.35 mm tele jack	
PUSH	x 1	-	Yes
PUSH+PEAR	x 1	6.35 mm tele jack	
PUSH+PULL	x 1	Pull cord	Yes
3PUSH	x 3		Yes
3PUSH+PEAR	x 3	6.35 mm tele jack	

ROOM Category products

ROOM	x 3	-	Yes
ROOM+PEAR	x 3	6.35 mm tele jack	

4 Hardware overview

4.1 Front overview

The front consist of a plastic cover with one or three LEDs and a combination of buttons and connectors (examples below).

PUSH+PEAR

ROOM+PEAR/3PUSH+PEAR


KEY

4.2 Back overview

On the back are the connectors and contacts.


Denomination

1	Place for 6.35 mm tele jack connector/Pull cord/Key cylinder
2	Jumper J2
3	Jumper J1
4	AAA-battery connectors
5	Screw connector
6	Product (data) label

4.3 The screw connector

The screw connector allows for connecting external power, inputs and outputs by cable to the WALL unit.

4.4 Power

The unit can be powered either by 10-28V_{DC} or 2 x AAA batteries.

4.5 Inputs IN1 and IN2

In the screw connector there are two wired inputs, IN1 (pin 6) and IN2 (pin 5). Both inputs have a common ground (pin 7) and must be connected to potential free contacts.

Input activation is ignored during the first 10 seconds after start up.

4.6 Outputs OUT₁ and OUT₂

The outputs are galvanically separated from the wired inputs and each other and works according to below.

4.6.1 ROOM and ROOM+PEAR

The relay outputs act according to the settings in the builtin actionplan or by a superior system (by radio).

The builtin actionplan activates the outputs as below.

Mode	OUT ₁	OUT ₂
Standby	-	-
Alarm	Active	-
Presence	-	Active
Assistance	Active	Active

4.6.2 All other products

The relay outputs act according to settings made by manual configuration or by a superior system (by radio).

4.7 Jumpers J1 and J2

As default the jumpers **J1** and **J2** has the functions as described below.

	Mounted	Not mounted
J1	Receiver on, 24V mode*	Receiver off, battery mode
J2	Inputs (IN1 and IN2) normally closed**	Inputs (IN 1 and IN2) normally open**

However, jumper **J2** can be overridden by using WALL family Programmer (except for ROOM and ROOM+PEAR).

4.8 Buttons and LEDs

The buttons are used to trig alarms. For ROOM and ROOM+PEAR it will also change the unit's mode.

The different modes are:

- Standby mode
- Alarm mode
- Presence mode
- Assistance mode

The LEDs indicate action and/or the unit's state and can either blink or light steady depending on how the units is configured and/or used.

* If Jumper J1 is mounted the unit must be externally powered.

** IN1 for ROOM and ROOM+PEAR is always normally open.

4.9 Ports 1 and 2

The tele jack connector provides two ports when a stereo tele plug is used or one port if a mono teleplug is used .

The ports are of type **Normally Open**.

Both ports have a common ground and must be connected to potential free contacts (see below).


If an inserted tele plug is removed, the unit sends a **Tamper alarm**. If the plug is inserted again the unit sends a **Tamper OK alarm**.

4.10 Mounting frames

There are 3 different mounting frames to fix the unit on the installation. The IP44 mounting frame enables IP44 class protection for an IP44 ready unit, see “3 Product variant denominations”.


4.10.1 1 mm mounting frame

The 1 mm mounting frame is included in units by default and is used when mounting WALL in an electric coupling box.


4.10.2 13 mm mounting frame (accessory)

The standard 13 mm mounting frame fits all units and is used when the device is mounted on the wall surface.

4.10.3 Mounting frame IP44 (accessory)

The IP44 mounting frame is designed to provide IP44 class protection to IP44 ready product variants, see “3 Product variant denominations”. The IP44 mounting frame is also 13 mm high.


It is important to secure that the surface of the wall where the mounting frame is placed is flat and smooth so the mounting frame closes tightly to the wall.

The three mounting frames contain a plastic flange (marked with a red circle in the picture above) and is designed to fit on a new WALL PCB.

To determine if the WALL PCB is new or old, take a look at the picture below. New PCBs have a hole where the red rectangle is, while old PCBs do not.

i The mounting frame **can** be fitted on an old PCB, but then the plastic flange must be cut off.

5 Manual configuration

(not valid for ROOM and ROOM+PEAR)

Any unit can be manually configured regarding adding transmitters and how the two relay outputs should act when the added transmitter is activated.

The procedure is in short:

- Determine what position to activate
- Add the radio transmitter
- Determine the action when the added transmitter is activated

Please read the following instruction carefully before attempting to manually configure a radio transmitter.

1. Power off the unit (remove both 24V_{DC} and batteries).
2. Remove **J2**.
3. Power on the unit and replace the J2 within 5 seconds from power on.
4. The LED blinks green 3 times.
5. After the initial “three green blink”, the unit now starts to blink red. The number of blinks indicates the current radio position, i.e. 3 red blinks equals radio position 3.

$$\overset{1}{*} \rightarrow \overset{2}{*} \rightarrow \overset{3}{*} \rightarrow \overset{4}{*} \rightarrow \overset{5}{*} \rightarrow \overset{6}{*} \rightarrow \overset{7}{*} \rightarrow \overset{8}{*} \rightarrow \overset{1}{*} \rightarrow \overset{2}{*}$$
6. To select radio position, press any button, port or input.
7. The unit confirms the selection by lighting green for 2 seconds.

8. Now activate the radio transmitter (e.g. press the red button on an SMILE).
9. The unit will confirm the successful reception of the radio ID code by lighting green for 2 seconds.
10. The unit now starts blinking green, allowing the setting the action for relay 1.

1 → 2 → 3 → 4 → 5 → 6 → 1 → 2
 * → * → * → * → * → * → * → *

11. The number of blinks indicates what action to set.

1 blink	Pulse
2 blinks	On
3 blinks	Off
4 blinks	Toggle
5 blinks	No action
6 blinks	Stop pulse

12. Press any port, button or input at the appropriate action blink.
13. The unit now starts blinking red, allowing the setting the action for relay 2 (see #10).

1 → 2 → 3 → 4 → 5 → 6 → 1 → 2
 * → * → * → * → * → * → * → *

14. Activate any port, button or input at the appropriate action blink.

15. If successful, the LED rapidly blinks green for 2 seconds.

 *At any stage during the manual configuration, if there is an unsuccessful configuration, the LED blinks rapidly for 2 seconds. When this happens the manual configuration is cancelled and it is necessary to restart the procedure.*

In the event of inactivity during manual configuration (e.g. a transmitter is not activated after #8 or a port/button is not pressed after #11) the unit waits for 60 seconds and then returns to normal mode. The manual configuration must then be restarted.

5.1 Remove a transmitter

To remove a transmitter, follow steps #1 through #7 above. After step #7, remove jumper **J2**. The device will confirm the deletion by rapidly blinking green for 2 seconds.

5.2 WALL family Programmer

All WALL units can be configured with a Neat Programming Unit (NPU) and the software WALL Family Programmer.

The NPU can be purchased from NEAT Electronics and the software can be downloaded from our homepage:
www.neat-group.com/downloads/documentation

5.3 More information

For more information, please refer to WALL Family - Technical Handbook, NE41 11014-02.

6 Technical data

Dimension (excl. wallframe)	86 x 86 x 14 mm
Weight (incl. batteries)	64 g
Supply voltage _{min}	10 V _{DC}
Supply power _{max}	28 V _{DC} , 100 VA
Battery type	GP Ultra Alkaline AAA or GP Super Alkaline AAA
Battery life*	5 years
Frequency _{RF} - EU	869.200 - 869.250 MHz (Social alarms) Category 1 869.400 - 869.450 MHz (Social alarms) Category 2
Frequency _{RF} - US	916.2, 916.4 MHz
Frequency _{RF} - Other**	866, 868, 906, 916, 921 MHz
Transmit power _{max} - EU	869.2 MHz - 10 mW (EIRP) 869.4 MHz - 500 mW (EIRP)
Transmit power _{max} - Other	According to local regulations
Relay output, switching _{max}	30 V _{DC} , 250 mA
Button push time _{min}	150 ms
Wired input open/close _{min}	150 ms
Jumper change delay	3 s
Tamper remove/insert delay	3 s

 *The frequencies for each market are programmed during production and not changeable by the user.*

* Radio test alarm + 4 alarms every day (when battery powered).

** According to local regulations for social alarms.

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Die NPU und die Software kann bei NEAT GmbH erworben werden.